

## Ingredients

Refined Sunflower Oil, AB-KOLICARE® Probiotic Strains\* (*Pediococcus pentosaceus* CECT 8330 and *Bifidobacterium longum* CECT 7894).

## Instructions

**SHAKE WELL BEFORE USE**



The recommended daily dose is 5 drops per day, before or after the first meal. It is recommended to administer the product directly to the baby's mouth.

### Alternative options include:

Add 5 drops into baby formula.

Add 5 drops to clean finger and allow baby to suck. Add 5 drops to Mothers nipple and allow baby to suck.

Quantity of combined strains per daily dose:  $1 \times 10^9$  CFU. For children, drops can be mixed with milk (maximum temperatures 37°C).

**Do not mix with boiling milk.**

**SAFE:** Bacteria originating from human breastmilk and intestinal microbiota of healthy babies. QPS and Self Affirmed GRAS strains.

## Features **100% NATURAL\***



100% Organic Strains\*



Suitable from Birth



Good for the Gut



Lactose Free



GMO Free



Gluten Free



For Babies, Infants & Children



Formula & Breastfed Babies



Free from Sugar, Alcohol & Colourant

**100% Organic Strains\***  
Suitable from Birth



**Clinically Proven\***

\*Clinical studies and certificates available on request.

# KoliCare

## For Babies, Infants and Children

**Clinically proven to reduce episodes of colic & babies daily crying time\***

Contains **Good Bacteria** for your baby's gut

## Storage

Store away from light in a cool, dry place not above 25°C.

Do not freeze. Once opened, consume within 1 month.

## Precautions

Food supplements cannot replace a varied and balanced diet or a healthy lifestyle. Do not exceed the recommended daily dose. Keep out of reach of young children.

**Best before** date see bottom of the box.

**Guaranteed Irish**

**MyPro KoliCare a ZEON Healthcare Ltd. brand.**

**Distributed by:**  
Pharmed Limited,  
Unit 3 Clonmore Business Park,  
Mullingar, Co. Westmeath,  
N91 DK59, Ireland.


**Licensed by:**  
AB-BIOTICS, S.A., Esade Creapolis,  
Av. de la Torre Blanca, 57 - Desp.  
3B11, 08172 St. Cugat del Valles  
(Barcelona), Spain.

Manufactured in the EU

8ml e

# HOW TO USE KoliCare


## Directions of use



1

### STEP 1


Remove the cap of the glass bottle (A) and replace it for the administration cap (B+C) screwing it well (not too hard!). You can store the product this way.



2

### STEP 2


Before administration, shake vigorously.



3

### STEP 3


Open the bottle by placing the fingers in the striated area of the cap (C) and twisting slightly, counter-clockwise.



4

### STEP 4


Apply the product by turning the bottle and pinching the plastic dropper (B) as many times as drops needed (daily dose: 5 drops).



5


### STEP 5

Close the bottle by adding the cap (C) and slightly twisting clockwise. Store until next administration.



6

**To clean** the plastic dropper (B), remove the cap (C) and unscrew the plastic dropper (B) by placing the fingers in the white area. You can wash it with water or boil it to sterilize it.

 5 Drops per day

 Easy to use

 Removable Dropper, Easy to Clean

 No Refrigeration Required